

***Concerned About The Serious
Consequences Of Peri-implantitis?***

*BIOMET 3i's Prospective, Multicenter, Randomized-Controlled Five-year Study Of
Hybrid OSSEOTITE® And Full OSSEOTITE Implants For The Incidence Of Peri-implantitis*

Peri-implantitis: Potential For Implant Failure


Peri-implantitis presents as a potentially serious clinical problem for patients and clinicians. It also impacts the viability of the dental implant as a treatment option for missing teeth. It can be a prominent cause for late implant failure leading to loss of the prosthesis.

What Is Peri-implantitis?


Peri-implantitis is a syndrome characterized by three clinical findings.

- **Severe mucosal inflammation (mucositis)**
- **Marked soft tissue Clinical Attachment Loss (CAL)**
- **Progressive crestal bone regression**


In order for a case to be declared as peri-implantitis, all three must be present with a primary microbial etiology.¹



Occlusal image of the implant (mandibular left side), reflecting suppuration and severe tissue inflammation due to microbial infection.



Periapical radiograph of the implant in the above clinical image. Significant bone loss is evident around the implant. The implant had +3 mobility and required removal and replacement.




Large peri-implantitis defect undergoing revision. Note the distance of the defect margin to the implant surface.


Peri-implantitis And Roughened Surfaces

Peri-implantitis is difficult to treat and may often lead to progressive bone loss and implant failure. Implants with roughened surfaces on coronal-implant collars may be perceived as having higher risks of peri-implantitis or, at the very least, other mucosal complications.

The incidence of peri-implantitis has been reported to be as high as 14.4%² and the prevalence as high as 28%.³ The risk of peri-implant disease had been thought to increase with greater implant surface roughness. Historically, very roughened implants (TPS and HA legacy coated implants) were reported as having improved initial integration success,⁴ but also associated with a higher proportion of late failures, some due to peri-implantitis.⁵



Patient case demonstrating peri-implantitis around Titanium Plasma Sprayed (TPS) Implants



SEM of TPS surface at 2000x magnification

Concerns about implant failure remain with roughened implants. Is this perception a reality for all roughened-surface implants?

One-Of-A-Kind Study

Addressing concerns about peri-implantitis and roughened-surface implants.

Historically, a machined surface implant has been recognized for its ability to be decontaminated, compared with roughened surfaces.^{6, 7, 8}


Acknowledging a clinical concern about the occurrence of peri-implantitis, BIOMET *3i* initially offered the OSSEOTITE® Implant with a hybrid surface design where the implant was machined from the abutment seating platform to the third thread and dual acid etched (DAE) to the apex.

The potential benefit of having the dual acid etched OSSEOTITE Surface complexity along the entire length of the implant was considered and developed. Yet, the question remained: How would the benefits of this dual acid etched surface play


against the possibility of increasing the incidence of peri-implantitis? This led to a specific effort to quantify the risk of adverse events for fully-etched implants vs. the hybrid surface design.

Recognizing the responsibility as an implant manufacturer to pursue evidence-based research, BIOMET *3i* sponsored a prospective, randomized-controlled clinical trial to determine if a difference exists in the incidence of peri-implantitis between hybrid and fully-etched implants.⁹

After a review of the published literature, this is the only prospective, randomized-controlled study with peri-implantitis as the primary outcome.



The control OSSEOTITE Implant is dual acid etched (DAE) from the apex to the third coronal thread. A machined surface continues to the seating platform.



The test OSSEOTITE Implant is dual acid etched (DAE) from the apex to the abutment seating platform. Both test and control implants are cpTi with straight walls, apical cutting features and an external hex connection.

Study Outcomes

A total of 304 implants were placed, supporting 127 prostheses, with a distribution of 139 control and 165 test implants in 112 patients.

Follow-up evaluations included:

- **Sulcus Bleeding Index Scores (SBI)**
- **Probing for suppuration**
- **Assessments for mobility**
- **Serial Periapical radiographs to identify radiolucencies and crestal bone levels**

Sulcus Bleeding Index

		Hybrid Control (%)	FOSS Test (%)
SBI Scores	0	83.5	84.3
	1	13.6	13.1
	2	2.6	2.4
	3	0.3	0.2

84% of all SBI scores were “0” (absence of bleeding); 13% of scores were “1” - isolated bleeding spot for both Full OSSEOTITE® (FOSS) and hybrid OSSEOTITE Implants.

Probing Depth Scores (Number Of Sites Probed)

		Control (N)	Test (N)
Probing Depths: Change From Baseline (mm)	0 ≤ 1	147	119
	1.1 ≤ 3	36	35
	3.1 ≤ 5	0	0
	> 5	0	0

No implant (test or control) showed changes in probing depths greater than 3.0mm.

No substantial differences in mucosal health outcomes between test and control groups were observed throughout the 5-year follow-up.

Only one observation of suppuration was recorded and it was for a control implant at the baseline evaluation. There was one diagnosis of peri-implantitis for a control implant 3.5 years after implant placement.

Radiographic analyses of crestal bone recession demonstrate that the mean change from baseline (provisionalization) is less for test implants in comparison to control implants (P<.0001).

Regressive Bone Remodeling



Full OSSEOTITE (FOSS) averaged less bone regression vs. hybrid OSSEOTITE Implants over the five year period of follow-up.

These findings are consistent with previous studies showing that the OSSEOTITE Implant Surface had no difference in soft tissue response when compared to a machined surface.^{10, 11}

For dental implants, a combination of optimal osseous fixation properties and a low risk for peri-implantitis are desired. The OSSEOTITE Surface has more than a decade of clinical use and evidenced-based research to support its efficacy. The results of this multicenter study show no increased risk in soft tissue complications or peri-implantitis for the studied FOSS Implants.


Hybrid OSSEOTITE® Vs. Full OSSEOTITE Surface Comparison

Full OSSEOTITE (FOSS) Implants have increased OSSEOTITE Surface area as compared to a hybrid OSSEOTITE Implant as shown below. With increased surface area, the FOSS Implant provides more contact for osteogenesis.

Scenarios in which FOSS Implants could be used:

- **Poor quality bone**
- **Short implant cases**
- **Sinus graft cases**
- **Immediate loading cases**
- **Platform switching cases**

Representative Differences In OSSEOTITE Surface Area



Implant	% More OSSEOTITE Surface Area
FOS410	32.1%
FOS510	51.7%
FOS685	60.9%

References

1. Mombelli A. Microbiology and antimicrobial therapy of peri-implantitis. *Periodontol 2000*. 2002;28:177-189.
2. Berglundh T, Persson L, Klinge B. A systematic review of the incidence of biological and technical complications in implant dentistry reported in prospective longitudinal studies of at least 5 years. *J Clin Periodontol*. 2002;29(SUPPL):197-212, 232-233.
3. Zitzmann NU, Berglundh T. Definition and prevalence of peri-implant diseases. *J Clin Periodontol*. 2008 Sep;35(8 Suppl):286-91.
4. Watson CJ, Tinsley D, Ogden AR, Russell JL, Mulay S, Davison EM. A 3 to 4 year study of single tooth hydroxylapatite coated endosseous dental implants. *Br Dent J*. 1999 Jul 24;187(2):90-4.
5. Rosenberg, et al. A comparison of characteristics of implant failure and survival in periodontally-compromised and periodontally-healthy patients: A clinical report. *Int J Oral Maxillofac Implants V19*. 2004:873-879
6. Dennison DK, Huerzeler MB, Quinones C, Caffesse RG. Contaminated implant surfaces: an in vitro comparison of implant surface coating and treatment modalities for decontamination. *J Periodontol*. 1994 Oct;65(10):942-8.
7. Teughels W, Van Assche N, Sliepen I, Quirynen M. Effect of material characteristics and/or surface topography on biofilm development. *Clin Oral Implants Res*. 2006 Oct;17 Suppl 2:68-81.
8. Burgers R, Gerlach T, Hahnel S, Schwarz F, Handel G, Gosau M. In vivo and in vitro biofilm formation on two different titanium implant surfaces. *Clin Oral Implants Res*. 2010;21:156-164.
9. Zetterqvist L, et al. A Prospective, Multicenter, Randomized-controlled, Five-year Study of Hybrid and Fully-etched Implants for the Incidence of Peri-implantitis (The FOSS Study) *J Periodontology*. 2009;Dec 23: (14pp.) [Epub ahead of print]
10. Abrahamsson I, et al. The mucosal attachment to titanium implants with different surface characteristics. *J Clin Periodontol*. 2002: 29: 448-455.
11. Zitzmann NU et al. Soft tissue reactions to plaque formations at implant abutments with different surface topography. *J Clin Periodontol*. 2002;29: 456-461.

Ordering Information

Certain® Internal Connection I M P L A N T S Y S T E M

Full OSSEOTITE® Parallel Walled Certain Implants

Length	3.25mm (D)	4.0mm (D)	5.0mm (D)	6.0mm (D)
8.5mm	IFOSM385	IFOS485	IFOS585	IFOS685
10.0mm	IFOSM310	IFOS410	IFOS510	IFOS610
11.5mm	IFOSM311	IFOS411	IFOS511	IFOS611
13.0mm	IFOSM313	IFOS413	IFOS513	IFOS613
15.0mm	IFOSM315	IFOS415	IFOS515	IFOS615

Full OSSEOTITE XP® Certain Implants

Length	4/5mm (P)	5/6mm (P)
8.5mm	IFOS4585	IFOS5685
10.0mm	IFOS4510	IFOS5610
11.5mm	IFOS4511	IFOS5611
13.0mm	IFOS4513	IFOS5613
15.0mm	IFOS4515	IFOS5615

Full OSSEOTITE Tapered Certain Implants

Length	3.25mm (D)	4.0mm (D)	5.0mm (D)	6.0mm (D)
8.5mm	IFNT3285	IFNT485	IFNT585	IFNT685
10.0mm	IFNT3210	IFNT410	IFNT510	IFNT610
11.5mm	IFNT3211	IFNT411	IFNT511	IFNT611
13.0mm	IFNT3213	IFNT413	IFNT513	IFNT613
15.0mm	IFNT3215	IFNT415	IFNT515	IFNT615

OSSEOTITE Certain PREVAIL® Implants

Length	3/4/3mm (P)	4/3mm (P)	4/5/4mm (P)	5/4mm (P)	5/6/5mm (P)
8.5mm	IIOS3485	IIOS4385	IIOS4585	IIOS5485	IIOS5685
10.0mm	IIOS3410	IIOS4310	IIOS4510	IIOS5410	IIOS5610
11.5mm	IIOS3411	IIOS4311	IIOS4511	IIOS5411	IIOS5611
13.0mm	IIOS3413	IIOS4313	IIOS4513	IIOS5413	IIOS5613
15.0mm	IIOS3415	IIOS4315	IIOS4515	IIOS5415	IIOS5615

External Hex Connection I M P L A N T S Y S T E M

Full OSSEOTITE Parallel Walled Implants

Length	3.25mm (D)	3.75mm (D)	4.0mm (D)	5.0mm (D)	6.0mm (D)
7.0mm	FOSM307	FOS307	FOS407	FOS507	FOS607
8.5mm	FOSM385	FOS385	FOS485	FOS585	FOS685
10.0mm	FOSM310	FOS310	FOS410	FOS510	FOS610
11.5mm	FOSM311	FOS311	FOS411	FOS511	FOS611
13.0mm	FOSM313	FOS313	FOS413	FOS513	FOS613
15.0mm	FOSM315	FOS315	FOS415	FOS515	FOS615

Full OSSEOTITE XP Implants

Length	3/4mm (P)	4/5mm (P)	5/6mm (P)
7.0mm	FOS3207	FOS4507	FOS5607
8.5mm	FOS3285	FOS4585	FOS5685
10.0mm	FOS3210	FOS4510	FOS5610
11.5mm	FOS3211	FOS4511	FOS5611
13.0mm	FOS3213	FOS4513	FOS5613
15.0mm	FOS3215	FOS4515	FOS5615

Full OSSEOTITE Tapered Implants

Length	3.25mm (D)	4.0mm (D)	5.0mm (D)	6.0mm (D)
8.5mm	FNT3285	FNT485	FNT585	FNT685
10.0mm	FNT3210	FNT410	FNT510	FNT610
11.5mm	FNT3211	FNT411	FNT511	FNT611
13.0mm	FNT3213	FNT413	FNT513	FNT613
15.0mm	FNT3215	FNT415	FNT515	FNT615



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside The U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com

To Receive Information About BIOMET 3i's Products, Services And Events
By Email, Please Visit Our Website At www.biomet3i.com



REV A 03/10

Certain, OSSEOTITE, OSSEOTITE XP and PREVAIL are registered trademarks of BIOMET 3i LLC. BIOMET is a registered trademark and BIOMET 3i and design are trademarks of BIOMET, Inc. ©2010 BIOMET 3i LLC. All rights reserved.

SUBSIDIARIES

AUSTRALIA
Phone: +61-2-9855-4444
Fax: +61-2-9888-9900

AUSTRIA
Phone: +43-(0)6235-200-45
Fax: +43-(0)6235-200-45-9

BELGIUM
Phone: +32-2-5410290
Fax: +32-2-5410291

BRAZIL
Phone: +55-11-5081-4405
Fax: +55-11-5081-7484

CANADA
Phone: +514-956-9843
Fax: +514-956-9844

FRANCE
Phone: +33-1-41054343
Fax: +33-1-41054340

GERMANY
Phone: +49-721-255177-10
Fax: +49-721-255177-73

IRELAND
Phone: +353-1-800-552-752
Fax: +44-1628-820182

JAPAN
Phone: +81-66868-3012
Fax: +81-66868-2444

KOREA
Phone: +82-2-5678-550
Fax: +82-2-5678-577

MEXICO
Phone: +52-55-2282-0120
Fax: +52-55-2282-0120 ext. 20

THE NETHERLANDS
Phone: +31-(0)78-629-2800
Fax: +31-(0)78-629-2801

NEW ZEALAND
Phone: +64-508-122-221
Fax: +64-508-133-331

NORDIC REGION
Phone: +46-40-17-6090
Fax: +46-40-17-6099

PORTUGAL
Phone: +351-21-000-1647
Fax: +351-21-000-1675

SPAIN
Phone: +34-93-470-59-50
Fax: +34-93-372-11-25

SWITZERLAND
Phone: +41-(0)44-200-76-00
Fax: +41-(0)44-200-76-01

U.K.
Phone: +44-800-652-1233
Fax: +44-1628-820182

DISTRIBUTORS

ARGENTINA
Dentalmax, SA
Phone: +541-1482-71001
Fax: +541-1482-67373

CHILE
Cybel, SA
Phone: +56-2-2321883
Fax: +56-2-2330176

CHINA
Atek Inc.
Phone: +86-21-6329-1265
Fax: +86-21-6329-1620

COLOMBIA
3i Colombia
Phone: +571-612-9362
Fax: +571-620-6412

COSTA RICA
Implantec S.A.
Phone: +506-234-9043
Fax: +506-224-7620

EL SALVADOR
Dentimer SA de CV
Phone: +503-263-6350
Fax: +503-263-6676

GREECE
Impladent Dental Implants, LLC
Phone: +30-2310-501-651
Fax: +30-2310-862-090

ISRAEL
H.A. Systems
Phone: +972-3-6138777
Fax: +972-3-6138778

HONG KONG
Ositek Inc., Ltd.
Phone: +852-8121-6601
Fax: +852-3747-3754

ITALY
Biomax, srl.
Phone: +39-0444-913410
Fax: +39-0444-913695

MIDDLE EAST
3i MENA s.a.l.
Middle East And North Africa
Phone: +961-1-694000
Fax: +961-1-694222

PARAGUAY
Andres H. Arce y Cia SRL
Phone: +595-21-208185
Fax: +595-21-496291

POLAND
Dental Depot Wasio
Phone: +48 71 335 70 71
Fax: +48 71 335 70 90

RUSSIA
Com-Dental
Phone: +7-495-580-3080
Fax: +7-495-580-3081

SINGAPORE
Asia Implant Support & Services
Phone: +65-6223-2229
Fax: +65-6220-3538

SOUTH AFRICA
Selective Surgical CC
Phone: +27-11-991-7007
Fax: +27-11-672-1391

TAIWAN
Kuo Hwa Dental Suppliers Co., Ltd.
Phone: +886-2-2226-1770
Fax: +886-2-2226-8747

THAILAND
3i (Thailand) Co., LTD.
Phone: +662-252-6685
Fax: +662-252-6686

TURKEY
Tamer Med. A.Ş.-Turkey
Phone: +90-212-465-3352
Fax: +90-212-465-3502

UKRAINE
Com-Dental
Phone: +38-067-7007667
Fax: +38-044-5017117

URUGUAY
Pro3i Implant S.R.L.
Phone: +598-2-4034163
Fax: +598-2-4034163

VIETNAM
D.O.E.
(Dentistry Of Excellence Vietnam)
Phone: +84 89 25 37 03
Fax: +84 89 25 37 00