

New as of:

08.2018



CEREC Optispray

Operating Instructions (valid for USA)

English



Table of contents

1	General data	3
1.1	Symbols used	3
1.2	Contact information	4
1.3	General conventions and structure of the document	5
1.3.1	Structure of the document.....	5
1.3.1.1	Identification of the danger levels	5
1.3.1.2	Formats and symbols used	5
2	Safety instructions	6
3	Product description	7
3.1	Composition	7
4	Intended use, indications and contraindications	8
4.1	Intended use	8
4.2	Indication.....	8
4.3	Contraindications	8
5	Application procedure	9

1 General data

1.1 Symbols used



CE mark

This product bears the CE mark in accordance with the provisions of the Council Directive 93/42/EEC of June 14, 1993 concerning medical devices (MDD).

USA: Rx only

For the USA only

CAUTION: Federal law (USA) restricts sale of this device to or on the order of a physician, dentist, or licensed practitioner.

REF ABC123

Legend

Article number

LOT ABC123

Batch number



Observe accompanying documents



The accompanying documents are available on the Dentsply Sirona homepage.

www.dentsplysirona.com/manuals

 YYYY-MM-DD

Date of manufacture: Year-Month-Day

 YYYY-MM-DD

Can be used until: Year-Month-Day

XXX 3

Filling volume of packaging

XXX ml e

Nominal filling quantity




Keep out of direct sunlight



Irritant



Easily or highly flammable



Temperature during storage and transport

1.2 Contact information

Customer service center

For technical questions, use the contact form on the internet at the following address:

<http://srvcontact.sirona.com>

Manufacturer's address



Sirona Dental Systems GmbH
Fabrikstrasse 31
64625 Bensheim
Germany

Tel.: +49 (0) 6251/16-0

Fax: +49 (0) 6251/16-2591

e-Mail: contact@dentsplysirona.com




www.dentsplysirona.com

1.3 General conventions and structure of the document

1.3.1 Structure of the document

1.3.1.1 Identification of the danger levels

To prevent personal injury and material damage, please observe the warning and safety information provided in these operating instructions. Such information is highlighted as follows:

 DANGER
An imminent danger that could result in serious bodily injury or death.
 WARNING
A possibly dangerous situation that could result in serious bodily injury or death.
 CAUTION
A possibly dangerous situation that could result in slight bodily injury.
NOTE
A possibly harmful situation which could lead to damage of the product or an object in its environment.
IMPORTANT
Application instructions and other important information.

Tip: Information for simplifying work.

1.3.1.2 Formats and symbols used

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"> ✓ Prerequisite 1. First action step 2. Second action step or ➤ Alternative action ↔ Result ➤ Individual action step 	Prompts you to do something.
See "Formats and symbols used [→ 5]"	Identifies a reference to another text passage and specifies its page number.
• List	Designates a list.
"Command / menu item"	Indicates commands / menu items or quotations.

2 Safety instructions

IMPORTANT

Shake before using!

Failure to observe the specified application restrictions and procedure may result in measuring errors.

CAUTION

No side effects are known to date. Allergic reactions may occur in individual cases. CEREC Optispray should not be used if it is suspected or known that the patient may be allergic to any of its ingredients.

CEREC Optispray must not be used on infants. Infants must not be present during a CEREC Optispray treatment.

3 Product description

CEREC Optispray enables easy and perfect optical impressions to be taken of CAD/CAM restorations. CEREC Optispray is designed for (chairside) intraoral application. Due to its optimized dispensing and a special nozzle, it always ensures homogeneous coating of the object to be recorded.

3.1 Composition

Pigment suspension in ethanol; fluorinated hydrocarbon.

Flavoring: Peppermint

4 Intended use, indications and contraindications

4.1 Intended use

CEREC Optispray is used as the coating for optical impressions with the CEREC system. It is used for the intraoral measurement of prepared and unprepared teeth and the surrounding area.

4.2 Indication

Intraoral application (prepared and unprepared teeth and the surrounding area)

4.3 Contraindications

- on uncleaned surfaces,
- on wet surfaces,
- on heavily bleeding wounds that cover a large area.

5 Application procedure

1. The surface to be coated must be thoroughly cleaned and dried.
2. Take the spray head with cannula/nozzle and slide it into the spray head until the colored marking on the cannula/nozzle disappears inside the spray head. Before every application, the secure fit of the cannula/nozzle must be checked by pulling it gently.
3. Shake the can well before use.
4. The cannula/nozzle can be rotated as required to enable optimal coating from all directions. When coating objects in the posterior tooth region, you can use the stabilizing tube to push the patient's lip to one side without deforming the cannula (the stabilizing tube is included only in the 200 ml package).
5. Use the spray jet to cover the exposure area.
6. The spray nozzle should be held approx. 10-15 mm away from the object.
7. Take an exposure with the 3D camera in the usual manner.
8. As soon as the optical 3D exposure has been taken, the surface must be cleaned with air/water spray.
9. The cannula/nozzle should be replaced after each use.

We reserve the right to make any alterations which may be required due to technical improvements.

© Sirona Dental Systems GmbH
D3344.201.04.11.02 08.2018

Sprache: englisch
Ä.-Nr.: 000 000

Printed in Germany

Sirona Dental Systems GmbH



Fabrikstr. 31
64625 Bensheim
Germany
www.dentsplysirona.com

Order No **66 87 698 D3344**