



The Protection You Need To Grow

OsseoGuard™ Resorbable Collagen Membrane



- Fibrillar Matrix Structure Provides Strength For Tacking Or Suturing If Desired
- Excellent Handling Characteristics When Hydrated Adaptable To Various Defects



Nurture The Environment For Predictable Results

Clinicians can rely on the OsseoGuard™ Resorbable Collagen Membrane to provide predictable outcomes and ease of use during Guided Bone Regeneration (GBR) procedures. Designed for optimal strength, resorption and handling, the OsseoGuard Membrane can help regenerate a variety of defect sites. A unique manufacturing process creates a longer resorption profile and the membrane's fibrillar matrix structure provides strength for tacking or suturing. The OsseoGuard Membrane is biocompatible for implantation and the pure bovine Type I collagen, sourced from closed herds (tracked from birth), increases security.

The following clinical indications describe where the OsseoGuard Membrane can be applied as a protective barrier to support the process of bone regeneration in patients:

- Localized Ridge Augmentation/Future Site Preparation
- Peri-Implant Bone Defects
- Extraction Sockets
- Bone Regeneration After Root Resection
- Sinus Window

Clinicians Who Are Seeing Growth



"BIOMET *3i*'s OsseoGuard Guided Bone Regeneration Barrier has helped me deliver bone growth with easier use than the other barriers currently available."

- Dr. Peter Shatz, United States



"The OsseoGuard Membrane reconstitutes quickly and handles well when used for socket grafting procedures or on the lateral wall of the maxillary sinus window."

- Dr. Lee Silverstein, United States



"The consistency and tear resistance are very good, even when wet."

- Dr. Roberto Cochetto, Italy



"The use of this kind of membrane retains the bone graft inside the sinus and the membrane's great stiffness protects the schneiderian membrane."

- Professor José Calvo, Spain



"The OsseoGuard Membrane has the best handling characteristics of any resorbable barrier I've used. Its stiffness makes it easy to insert and it remains in place with flap closure even without tacking."


- Dr. Paul Ricchetti, United States



OsseoGuard ™: A Case For Growing Bone



1. Pre-op radiograph shows a failing four-unit bridge and ridge defect.



2. Flap reflection reveals large ridge defects at #5 and 6 edentulous sites.



 Defects are grafted with freeze-dried mineralized bone allograft and covered with an OsseoGuard Membrane. Primary closure is achieved.



4. Healing after six weeks. Note excellent soft tissue response.



5. Re-entry at five months confirms significant bone regeneration.




6. Post-op radiograph at five months demonstrates a regenerated ridge.

Case Study Provided By: Dr. Michael Sonick

Verified Lot-to-Lot Strength

In order to evaluate the strength of the OsseoGuard Membrane, a suture pull-out test is conducted by the manufacturer (Collagen Matrix, Inc.) on every lot of membranes produced.

A 3.0 suture is passed through the membrane about 3mm from the edge of the membrane. A knot is tied in the suture, leaving a loop to hook the suture onto a force gauge. The other end of the membrane is secured in a clamp. The suture is pulled at a rate of one inch per minute until the suture pulls out. The average force required to pull a suture out of the OsseoGuard Membrane from ten recent suture pullout tests is: .286kg +/- 0.090kg.



Ordering Information

OsseoGuard™ Resorbable Collagen Membrane

Size: 15mm x 20mm

> Size: 20mm x 30mm

Catalog Number OG1520

Catalog Number OG2030



Catalog Number OG3040



The OsseoGuard Membrane Is Available In Three Convenient Sizes, Is Double Sterile Packaged And Has A Shelf-Life Of Three Years

For More Information On How To Grow Your Business With The OsseoGuard Membrane, Call Your BIOMET 3i Representative Today Distributed By:



Global Headquarters 4555 Riverside Drive Palm Beach Gardens, FL 33410 1-800-342-5454 Outside The U.S.: +1-561-776-6700 Fax: +1-561-776-1272

Sign Up For BIOMET 3i's Electronic Newsletter "BIOMET 3innovations." Simply Go Online To www.biomet3i.com

Manufactured By Collagen Matrix, Inc., Franklin Lakes, NJ

OsseoGuard is a trademark of BIOMET 31, Inc. BIOMET is a registered trademark and BIOMET 31 and design are trademarks of BIOMET, Inc. ©2008 BIOMET 3i, Inc. All rights reserved



REV A 01/08

SUBSIDIARIES

AUSTRALIA

www.biomet3i.com

Phone: +61-2-9855-4444 Fax: +61-2-9888-9900

AUSTRIA

Phone: +43-(0)6235-200-45 Fax: +43-(0)6235-200-45-9

BELGIUM

Phone: +32-2-5410290 Fax: +32-2-5410291

Phone: +55-11-5081-4405 Fax: +55-11-5081-7484

CANADA

FRANCE

Phone: +514-956-9843 Fax: +514-956-9844

Phone: +33-1-41054343 Fax: +33-1-41054340

GERMANY Phone: +49-721-255177-10 Fax: +49-721-255177-73

IRELAND

Phone: +35-31-477-3925 Fax: +35-31-402-9590

MEXICO

Phone: +52-55-5679-1619 Fax: +52-55-5684-8098

THE NETHERLANDS

Phone: +31-(0)78-629-2800 Fax: +31-(0)78-629-2801

NEW ZEALAND

Phone: +64-508-122-221 Fax: +64-508-133-331

NORDIC REGION

Phone: +46-40-17-6090 Fax: +46-40-17-6099

PORTUGAL

Phone: +351-21-000-1645 Fax: +351-21-000-1675

Phone: +34-93-470-59-50 Fax: +34-93-372-11-25

SWITZERI AND

Phone: +41-44-380-46-46 Fax: +41-44-383-46 55

Phone: +44-1628-829314 Fax: +44-1628-820182

DISTRIBUTORS

ARGENTINA

Dentalmax, SA Phone: +541-1482-71001 Fax: +541-1482-67373

CHILE

Cybel, SA Phone: +56-2-2321883 Fax: +56-2-2330176

CHINA

Atek Inc. Phone: +86-21-6329-1265 Fax: +86-21-6329-1620

COLOMBIA

3i Colombia Phone: +571-612-9362 Fax: +571-620-6412

COSTA RICA

Phone: +506-2-256411 Fax: +506-2-247620

EL SALVADOR

Dentimerc SA de CV Phone: +503-263-6350 Fax: +503-263-6676

GREECE

Kostas Kornisorlis and Co. Phone: +30-2310-501-651 Fax: +30-2310-522-417

HONG KONG

Ositek Inc., Ltd. Phone: +852-8121-6601 Fax: +852-3747-3754

ISRAEL

H.A. Systems Phone: +972-3-6138777 Fax: +972-3-6138778

ITALY

Biomax, srl. Phone: +39-0444-913410 Fax: +39-0444-913695

Implant Innovations Japan Phone: +81-66-868-3012 Fax: +81-66-868-2444

KOREA

Jungsan Biomed Corp. Phone: +82-2-516-1808 Fax: +82-2-514-9434

LEBANON

Middle East And North Africa Phone: +961-1-694000 Fax: +961-1-694222

PARAGUAY

Andres H. Arce v Cia SRL Phone: +595-21-208185 Fax: +595-21-496291

POLAND

Dental Depot Phone: +48-71-341-3091 Fax: +48-71-343-6560

RUSSIA

Com-Dental Phone: +7-495-797-6686 Fax: +7-499-242-9567

SINGAPORE

Asia Implant Support & Services Phone: +65-6223-2229 Fax: +65-6220-3538

SOUTH AFRICA

Selective Surgical CC Phone: +27-11-991-7007 Fax: +27-11-672-1391

Kuo Hwa Dental Suppliers Co., Ltd. Phone: +886-2-2226-1770 Fax: +886-2-2226-8747

THAILAND

3i (Thailand) Co., LTD. Phone: +662-252-6685 Fax: +662-252-6686

UKRAINE

Com-Dental Phone: +38-067-7007667 Fax: +38-044-5017117

URUGUAY

Pro3implant S.R.L. Phone: +598-2-4034163 Fax: +598-2-4034163