

Anything Is Possible Providing You Work As A Team

The Provide® Abutment

- Providing Superior Surgical And Restorative Flexibility To Better Serve Patients
- A Strategic Tool For Practice Growth
- Snap-On Impression Components For Ease Of Use
- Certain® QuickSeat® Connection Produces
 An Audible And Tactile Click Confirming Proper Seating
- Kit Packaging For Streamlined Case Processing



Optimal Patient Care Starts With A Strong Product

With the increasing demands for dental implant therapy, the Provide® Abutment offers both surgical and restorative clinicians more options and greater flexibility to meet these demands and better serve patients. It allows them to partner together for a simpler restorative and surgical experience.

The implant is placed at bone level and not transgingivally, therefore the final crown margin is part of the abutment, not part of the implant. Provide Abutments give clinicians four collar height options for unmatched flexibility. This means the final crown margin is determined by restorative indication, not surgical consequence.

Providing Surgical Flexibility

- No Need To Determine Final Crown Margin At The Time Of Surgery
- Multiple Collar Heights (1, 2, 3 and 4mm)
- No Need To Prepare Implant Collar Due To Tissue Recession
- Use In A One Or Two-Stage Surgical Protocol
- Opportunity For Platform Switching™ To Potentially Preserve Crestal Bone




Achieve Your Practice Growth Goals Faster

A change of strategy could put you in the winner's circle. When clinicians want to reach their practice growth goals, using the Provide Abutment can be ideal. The simple restorative protocol is often easier than traditional crown and bridge procedures, which may allow for more implant referrals and greater adoption of the implant option.

Because restorative results and component costs are more predictable, you may achieve your practice growth goals more consistently. The Provide Abutment simplifies implant dentistry and may allow for more predictable growth of that segment of your practice.

Start To Finish, BIOMET 3i Helps Restorations Stay On Track

Provide® Component Selection Flowchart



^{*} b: collar height selection depends on tissue depth



A Winning Product Portfolio For Simpler Restorations

Provide® Abutments

Platform		4mm Pos	st Height		5.5mm Post Height					
Collar Heig	hts > 1mm	2mm	3mm	4mm	1mm	2mm	3mm	4mm		
4.1mm	IPA4140		IPA4340	IPA4440	IPA4155	IPA4255	IPA4355	IPA4455		
5mm	IPA5140	IPA5240	IPA5340	IPA5440	IPA5155	IPA5255	IPA5355	IPA5455		
6mm	IPA6140	IPA6240	IPA6340	IPA6440	IPA6155	IPA6255	IPA6355	IPA6455		

İ	Impression Copings Non-Prepared Post Profiles			<i>Prepare</i> Profi					ared Post les	Prepared Post Profiles			
	Post Height	4.8mm	6.5mm	4.8mm	6.5mm	311	Post Height	4.8mm	6.5mm	4.8mm	6.5mm		
	4mm	PIC484	PIC654	PIC484H	PIC654H		4mm	PAA484	PAA654	PAAP48	PAAP6		
	5.5mm	PIC485	PIC655	PIC485H	PIC655H		5.5mm	PAA485	PAA655	PAAP48	PAAP65		
	Protection Caps						Waxing Sleeves						
	Profiles				31		Profiles						
	Post Height	4.8mm	6.5mm				Post Height	4.8mm	6.5mm				
	4mm	PPC484	PPC654				4mm	PWS48S	PWS65S	Single Units			
	5.5mm	PPC485	PPC655				5.5mm	PWS48M	PWS65M	Multiple Un	its		
	Temporary Cylinders Profiles						Certain® G		Screw JNIHG				

Provide Kits

4mm,5.5mm PUA48S PUA65S 4mm,5.5mm PUA48M PUA65M

4.1mm Platform			5mm Platform				6mm Platfo		
Provide Abutment	Abutment Placement Kit*	Abutment Restorative Kit**	Provide Abutment	Abutment Placement Kit*	Abutment Restorative Kit**		Provide Abutment	Abutment Placement Kit*	Abutment Restorative Kit**
IPA4140	PAK4140	PRK484	IPA5140	PAK5140	PRK654		IPA6140	PAK6140	PRK654
IPA4240	PAK4240	PRK484	IPA5240	PAK5240	PRK654		IPA6240	PAK6240	PRK654
IPA4340	PAK4340	PRK484	IPA5340	PAK5340	PRK654		IPA6340	PAK6340	PRK654
IPA4440	PAK4440	PRK484	IPA5440	PAK5440	PRK654		IPA6440	PAK6440	PRK654
IPA4155	PAK4155	PRK485	IPA5155	PAK5155	PRK655		IPA6155	PAK6155	PRK655
IPA4255	PAK4255	PRK485	IPA5255	PAK5255	PRK655		IPA6255	PAK6255	PRK655
IPA4355	PAK4355	PRK485	IPA5355	PAK5355	PRK655		IPA6355	PAK6355	PRK655
IPA4455	PAK4455	PRK485	IPA5455	PAK5455	PRK655		IPA6455	PAK6455	PRK655

Single Units

Multiple Units

Please Note: Provide Impression and Restorative Components are not compatible with the Straumann-ITI® Implant System.

^{*} Abutment Placement Kit includes the Provide Abutment, Protection Cap and Abutment Screw.

** Abutment Restorative Kit includes the Provide Impression Coping (non-prepared), Provide Impression Coping w/Holes (Prepared), Provide Waxing Sleeve (single unit and multi-unit) and Provide Abutment Analog (non-prepared) - The Provide Abutment Analog w/Pin (prepared) is NOT included.



Easy-To-Use Components Get You To The Finish Line









Team Up With Your BIOMET 3i Sales Representative Today To Provide A Winning Combination For Your Patients And Your Practice.



Global Headquarters 4555 Riverside Drive Palm Beach Gardens, FL 33410 1-800-342-5454 Outside The U.S.: +1-561-776-6700 Fax: +1-561-776-1272 www.biomet3i.com

Sign Up For BIOMET 3i's Electronic Newsletter "3innovations." Simply Go Online To www.3i-online.com/signup



Rev C 02/07

Certain, Gold-Tite, Provide and QuickSeat are registered trademarks and Platform Switching is a trademark of Implant Innovations, Inc. BIOMET, 31 and design are trademarks and BIOMET is a registered trademark of BIOMET, Inc. LOCATOR is a registered trademark of Zest Anchors, Inc. ©2007 BIOMET, Inc. All rights reserved.

SUBSIDIARIES

Phone: +61-2-9855-4444 Fax: +61-2-9888-9900

Phone: +32-2-5410290

Fax: +32-2-5410291

BRAZIL

Phone: +55-11-5081-4405 Fax: +55-11-5081-7484

CANADA

Phone: +514-956-9843 Fax: +514-956-9844

FRANCE

Phone: +33-1-41054343 Fax: +33-1-41054340

Phone: +49-721-255177-10 Fax: +49-721-255177-73

IRELAND

Phone: +35-31-477-3925 Fax: +35-31-402-9590

MEXICO

Phone: +52-55-5679-1619 Fax: +52-55-5684-8098

THE NETHERLANDS

Phone: +31-(0)78-629-2800 Fax: +31-(0)78-629-2801

NEW ZEALAND

Phone: 0508-122-221 Fax: 0508-133-331

NORDIC REGION

Phone: +46-40-17-6090 Fax: +46-40-17-6099

PORTUGAL

Phone: +351-21-000-1645 Fax: +351-21-000-1675

Phone: +34-93-470-59-50 Fax: +34-93-372-11-25

SWITZERLAND Phone: +41-1-3804646 Fax: +41-1-3834655

Phone: +44-1628-829314 Fax: +44-1628-820182

DISTRIBUTORS

Dentalmax, SA Phone: +541-1482-71001 Fax: +541-1482-67373

AUSTRIA

Wie**l**adent Phone: +43-7672-93901 Fax: +43-7672-93903

CHILE Cybel, SA Phone: +56-2-2321883 Fax: +56-2-2330176

Atek Inc. Phone: +86-21-6329-1265 Fax: +86-21-6329-1620

COLOMBIA 3i Colombia

Phone: +571-612-9362 Fax: +571-620-6412

COSTA RICA

Implantec S.A. Phone: +506-2-256411 Fax: +506-2-247620

FL SALVADOR

Dentimerc SA de CV Phone: +503-263-6350 Fax: +503-263-6676

GREECE

Kostas Kornisorlis and Co Phone: +302310-269-079 Fax: +302310-555-573

ISRAEL

H.A. Systems Phone: +972-3-6138777 Fax: +972-3-6138778

ITALY

Biomax, srl. Phone: +39-0444-913410 Fax: +39-0444-913695

Implant Innovations Japan Phone: +81-66-868-3012 Fax: +81-66-868-2444

KOREA

Jungsan Biomed Corp. Phone: +82-2-516-1808 Fax: +82-2-514-9434

LEBANON

Middle East And North Africa Phone: +961-1-694000 Fax: +961-1-694222

PANAMA

Odontomedica, S.A. Phone: +507-2-239622 Fax: +507-2-239621

PARAGUAY

Andres H. Arce y Cia SRL Phone: +595-21-208185 Fax: +595-21-496291

POLAND

Dental Depot Phone: +48-71-341-3091 Fax: +48-71-343-6560

RUSSIA

Com-dental Phone: +7495-797-6686 Fax: +7495-242-9567

SINGAPORE

Asia Implant Support & Services Phone: +65-6223-2229 Fax: +65-6220-3538

TAIWAN

Kuo Hwa Dental Suppliers Co., Ltd. Phone: +886-2-2226-1770 Fax: +886-2-2226-8747

THAILAND

3i (Thailand) Co., LTD. Phone: +662-252-6685 Fax: +662-252-6686

UKRAINE

Olimpex, Ltd. Phone: +38-044-209-5713 Fax: +38-044-482-3678

URUGUAY

Pro3implant S.R.L. Phone: +598-2-4034163 Fax: +598-2-4034163